



Carefully read the centrifuge instruction manual and the rotor instruction manual before using the rotor.

P32ST2 Swinging bucket rotor

Max. Speed	32,000 rpm
Max. R.C.F.	193,000 xg
Rotor nominal capacity	16 mL × 6 = 96 mL
K factor	216
Weight	7.2 kg



(1) Summary of tubes

Tubes & Bottles					Adapter & Bucket			Cap	Maximum Speed (rpm)	Maximum R.C.F. (xg)
Nominal Vol. (ml)	Part No.	Name	Size (Φ x L cm)	Net Vol. (mL)	Part No. Name	Size (Φ x L cm)	Material	Part No. Name		
16	S303277A (50pcs/set)	16PET tube	16 x 102	17.0	---	---	---	---	32,000	193,000
16	332563A (50pcs/set)	16PA tube	16 x 102	15.0	---	---	---	---	32,000	193,000
3	336708A (50pcs/set)	3PA tube	11 x 47	2.8	336701A (6pcs/set) 3S16 adapter		POM	---	24,500	77,500
4	329561A (50pcs/set)	4PA tube	11 x 60	3.7	336702A (6pcs/set) 4S16 adapter		POM	---	24,800	86,800
5	332245A (50pcs/set)	5PA tube	13 x 52	4.9	336703A (6pcs/set) 5S16 adapter		POM	---	24,800	80,900
3	S301599A (50pcs/set)	3PC tube	13 x 51	3.4	336703A (6pcs/set) 5S16 adapter		POM	---	24,800	80,900
6.5	329445A (50pcs/set)	6.5PA tube	13 x 64	5.8	336704A (6pcs/set) 6.5S16 adapter		POM	---	24,800	89,200
6.5	326415A (50pcs/set)	6.5PC tube	13 x 64	6.5	336704A (6pcs/set) 6.5S16 adapter		POM	---	24,800	89,200
4.7	S303141A (50pcs/set)	4.7PC tick-walled tube	13 x 64	4.5	336704A (6pcs/set) 6.5S16 adapter		POM	---	24,800	89,200

(2) Standard Accessories

Part No.	Name	Qty.
332563A	16PA tube (50pcs/set)	2 (100)
84521013 ^{*1}	NBR O-Ring 1AP20	18
452701	Tweezers	1
209642A	Bucket Stand	1
216715	Rotor Stand	1
84810601	Lubricant for Screw	1
483719	Silicone Grease (Vacuum Grease)	1
S205074	Tool Box (C)	1

*1 NBR O-Ring 1AP20 (P/N 84521013) is sold in 5pcs/set at the additional purchase.

(3) Optional Accessories

Part No.	Name	Qty.
84521013	NBR O-Ring 1AP20	5pcs/set